

Smart & simple temperature monitoring



Temperature



Humidity

This is how simple the TempBeacon data logger works:

1. Place TempBeacon in the transport box near the medication.



2. Download Tec4App
(Playstore / Appstore).

3. Scan QR code with Tec4App and read
TempBeacon data automatically.



4. Check the temperature in the app

TempBeacon

Scan code for
data read-out.

Tec4med
Smart technology inside.
www.tec4med.com

Safely refrigerate medications throughout your journey



Cold packs



Carrying bag



Temperature
Data logger



Designed for
hand luggage

How to prepare the Cube:

1. Remove the 3 cold packs and place them in the freezer at approx. -18°C for at least 8 hours.

2. Defrost the cold packs at room temperature [approx. 30 min.]. The indicator changes color when the packs are ready.

3. Place cold packs, medicines and TempBeacon in the Cube.



4. Close the Cube lid and the carrying bag. Do not open more often than necessary!

5. Monitor the temperature with the TempBeacon [instructions on the back].

Smarte & einfache Temperaturüberwachung



Temperature



Humidity

So einfach funktioniert der TempBeacon Datenlogger:

1. TempBeacon in die Transportbox nahe dem Medikament legen.



2. Tec4App downloaden (Playstore / Appstore).

3. QR-Code mit Tec4App scannen und TempBeacon Daten automatisch auslesen.



4. In der App die Temperatur überprüfen

TempBeacon

Scan code for data read-out.

Tec4med

Smart technology inside.
www.tec4med.com

Medikamente sicher kühlen während Ihrer gesamten Reise



Cold packs



Carrying bag



Temperature
Data logger



Designed for
hand luggage

So bereiten Sie den Cube vor:

1. Die 3 Kühlakkus herausnehmen und in der Gefriertruhe bei ca. -18°C für mindestens 8 Stunden herunterkühlen.

2. Kühlakkus bei Raumtemperatur antauen (ca. 30 min.). Der Indikator verfärbt sich sobald die Akkus bereit sind.

3. Kühlakkus, Medikamente und TempBeacon in den Cube legen.



4. Cube Deckel und Tragetasche schließen. Nicht häufiger als notwendig öffnen!

5. Mit dem TempBeacon die Temperatur überwachen (Anleitung auf der Rückseite).